

# Máscara protectora antibacterial reutilizable

## Antibacterial protection re-usable mask

UNE 0065 ( Ensayos según UNE 14683:2019 +AC:2019 )



**Respirabilidad Breathability**  
<40 Pa/cm<sup>2</sup>



**Eficacia Filtración Bacteriana**  
(BFE) > 95,4%



**Material hidrófugo y filtrante  
tricapa**  
*Three layer water-repellent and  
filtering material*



**Esterilización Rayos Gamma**  
*Gamma Ray Sterilization treatment*



**Tratamiento Antibacteriano Sanitized T99-19**  
*Antibacterial Treatment Sanitized T99-19*



**PreformaHealth**  
by Preformatex

**FABRICANTES DESDE 2002**

**MADE IN SPAIN. HECHO EN ESPAÑA**

### Instrucciones de Uso / *Directions For Use:*

**Niños con supervisión de un adulto / *Always adult supervision for children***



**No usar secadora**  
*Not use dryer*



**Lavar a 60°** *Washed at 60 °*



**No planchar**  
*Not ironing*



**Eliminación bacterias. Activación  
protección antibacteriana. 10/15  
segundos**

*Bacteria elimination .Antibacterial  
protection reactivation.10/15 seconds*



**Sólo uso personal**  
*Only for personal use*



[www.preformahealth.es](http://www.preformahealth.es)

**OEKO-TEX®**  
CONFIDENCE IN TEXTILES  
**STANDARD 100**  
2017OK0373 AITEX  
Tested for harmful substances.  
www.oeko-tex.com/standard100